



**NEW**




**K-Shield® Advantage Phlebotomy Sets**

- Low-cost design incorporates pre-attached needles for accessing all vacuum bottles
- Fast and simple one-step activation with an “audible” click
- Pre-attached 16G venipuncture needle
- 16G needle with backeye optimizes blood flow
- Needle with wing for improved needle control
- Roller clamp to control blood flow
- 36-inch tubing length
- Latex-free



**Standard Phlebotomy Sets**

- Sets are available with and without pre-attached 16G venipuncture needle
- No safety device



Code #	Description (25 per box / 2 boxes per case)	Length
BO-1000	Phlebotomy set with 16G vacuum bottle needle	36 inches
BO-16S1	Phlebotomy set with 16G vacuum bottle needle and pre-attached 16G venipuncture needle with backeye	36 inches
BA-16S1	K-Shield Advantage Phlebotomy set with 16G vacuum bottle needle and pre-attached 16G venipuncture needle with backeye and safety device	36 inches