



K-Shield® Port Access Infusion Sets

CE 0123

- **Low profile design** – streamlined design helps make dressing simple and efficient and improves patient comfort during long infusion times
- **Simple safety device activation** – clinician-friendly and easy to use; audible click lets you know the needle is locked and enclosed during activation
- **Non-coring siliconized needle** – minimizes irritation while the port is being accessed or when the needle is removed
- **Spacers on skin plate** – allows air flow and minimizes moisture between patient’s skin and underneath skin plate
- **Non-DEHP fluid path** – offers a safe, DEHP-free fluid path

**Non-DEHP
Fluid Path**

Microbore 8” Tubing (20 per box / 5 boxes per case)

Code #	Needle Gauge	Needle Length	Pinch Clamp Color
PS03090	19	3/4”	Brown
PS03091	19	1”	Brown
PS03000	20	3/4”	Yellow
PS03001	20	1”	Yellow
PS03020	22	3/4”	Black
PS03021	22	1”	Black

Microbore 11” Tubing with Needle-Free Injection Site (20 per box / 5 boxes per case)

Code#	Needle Gauge	Needle Length	Pinch Clamp Color
PS01N90	19	3/4”	Brown
PS01N91	19	1”	Brown
PS01N00	20	3/4”	Yellow
PS01N01	20	1”	Yellow
PS01N20	22	3/4”	Black
PS01N21	22	1”	Black

Standard 8” Tubing (20 per box / 5 boxes per case)

Code #	Needle Gauge	Needle Length	Pinch Clamp Color
PS04091	19	1”	Brown
PS04000	20	3/4”	Yellow
PS04001	20	1”	Yellow



Spacers Minimize Moisture

Siliconized Needle