




Standard Port Access Infusion Sets € 0123

- Low-cost port access sets available in a variety of models
- Non-DEHP fluid path
- Non-coring needle
- Latex-free
- Needle hub and wing design offers secure and accurate needle placement



Non-DEHP Fluid Path



Microbore 8" Tubing (25 per box / 4 boxes per case)

Code #	Needle Gauge	Needle Length	Pinch Clamp Color
PI03090	19	3/4"	Brown
PI03091	19	1"	Brown
PI03000	20	3/4"	Yellow
PI03001	20	1"	Yellow
PI03020	22	3/4"	Black
PI03021	22	1"	Black

Microbore 11" Tubing with Y-injection Site (25 per box / 4 boxes per case)

Code#	Needle Gauge	Needle Length	Pinch Clamp Color
PI01Y90	19	3/4"	Brown
PI01Y91	19	1"	Brown
PI01Y00	20	3/4"	Yellow
PI01Y01	20	1"	Yellow
PI01Y03	20	1 1/2"	Yellow
PI01Y20	22	3/4"	Black
PI01Y21	22	1"	Black

Standard 8" Tubing (25 per box / 4 boxes per case)

Code #	Needle Gauge	Needle Length	Pinch Clamp Color
PI04091	19	1"	Brown
PI04000	20	3/4"	Yellow
PI04001	20	1"	Yellow