




PORT ACCESS INFUSION SETS

with High Pressure Tubing





K-Shield® Port Access Infusion Sets with High Pressure Tubing



- **Simple safety device activation**
- **Needle allows for higher pressure** and flow rates for power injections ports
- **Enables contrast agents** to be power-injected
- **Needles provide the flexibility of use** for any therapy need, including power injection of contrast media
- **Power sticker notifies CT/MRI tech** of power injection capability and CT/MRI parameters
- **19 gauge and 20 gauge** needles withstand up to 5 mL/s power injection of contrast media
- **22 gauge needle** withstands up to 2 mL/s power injection of contrast media
- **Non-DEHP fluid path** – offers a safe, DEHP-free fluid path
- **Latex-free**
- **Non-coring siliconized needle**

Non-DEHP Fluid Path

Microbore 7" Tubing (20 per box / 5 boxes per case)


Code #	Needle Gauge	Needle Length	Pinch Clamp Color
PH03090	19	3/4"	Brown
PH03091	19	1"	Brown
PH03000	20	3/4"	Yellow
PH03001	20	1"	Yellow
PH03020	22	3/4"	Black
PH03021	22	1"	Black

Microbore 8" Tubing with Open-End Cap Injection Site (20 per box / 5 boxes per case)

Code#	Needle Gauge	Needle Length	Pinch Clamp Color
PH01X90	19	3/4"	Brown
PH01X91	19	1"	Brown
PH01X00	20	3/4"	Yellow
PH01X01	20	1"	Yellow
PH01X20	22	3/4"	Black
PH01X21	22	1"	Black



Before activation



After activation

KAWASUMI
 KAWASUMI LABORATORIES AMERICA, INC.
www.kawasumiamerica.com

1-800-KAWASUMI (529-2786)
 P.O. Box 24355 • Tampa, FL 33623-4355
 Phone: 813-630-5554 • Fax: 813-630-5033
 e-mail: orders@kawasumiamerica.com

KAWASUMI LABORATORIES, INC.
 Shinagawa Intercity Tower B, 2-15-2,
 Konan, Minato-ku, Tokyo 108-6109, Japan