



K-Shield® Advantage Winged Blood Collection Sets


- Provides a safe collection environment for both clinician and patient
- Integrated safety device eliminates extra steps for disposal
- Fast and simple one-step activation with an “audible” click
- Needle is covered during removal, ensuring confidence in safety

K-Shield® Advantage Winged Blood Collection Sets with Multiple Sample Luer Adapter


(50 per box / 4 boxes per case)

Code#	Needle Gauge	Wing Color	Needle Length	Tubing Length
DBMA-21G	21g	Green	3/4"	12"
DBMA-23G	23g	Light Blue	3/4"	12"
DBMA-25G	25g	Dark Blue	3/4"	12"

K-Shield Advantage with Activated Safety Device




Raised Tab



Flexible Wings with Textured Surface and Color Coding for Easy Identification

K-Shield Advantage Safety Device

Female Luer Lock Connector



Easy-Grip Hub

Self-Sealing Synthetic Needle Sheath

Superior Needle Sharpness

Concave and Convex Dimples

Winged Blood Collection Set

Kink Resistant Tubing

