


Standard Winged Blood Collection Sets with Multiple Sample Luer Adapter

CE 0123

- Low cost winged collection set with multiple sample luer adapter
- Self-sealing synthetic needle sheath
- Flexible wings with textured surface for secure grip
- Easy-grip hub fits tube holders with anti-stick devices
- Sharpest needle point*



| Code # | Description (50 per box / 4 boxes per case) | Needle Gauge | Wing Color | Needle Length | Tubing Length |
|----------|---|--------------|------------|---------------|---------------|
| DBM1-21G | Winged collection set with multiple sample luer adapter | 21g | Green | 3/4" | 12" |
| DBM1-23G | Winged collection set with multiple sample luer adapter | 23g | Light Blue | 3/4" | 12" |
| DBM1-25G | Winged collection set with multiple sample luer adapter | 25g | Dark Blue | 3/4" | 12" |

*NAMSA needle sharpness report available upon request.



Intermittent Injection Site



Intermittent Injection Site

- Low cost injection site
- Male luer connector
- Allows clinician to perform multiple injections when attached to a winged infusion set
- Latex-free

Intermittent Injection Site

(50 per box / 4 boxes per case)

| Code# | Description |
|--------|---|
| INT-01 | Injection site with male luer connector |